

**NEMIS Special Webinar**


**27. Nov. 2024**




# Sampling of large surfaces with the N-Light™ MaxiSampler

**Julian Ihssen, NEMIS Technologies AG**

# Sampling: A crucial part of environmental monitoring in food factories



# Microorganisms on surfaces: an incredibly complex topic



- Strength of adhesion influenced by properties of surface [material, morphology] and cells [outer polysaccharides, protein fimbriae]
- Drying will increase adhesion
- Cells can be in resting state or injured [hard to recover]
- Any additional matrix [food residues] will affect survival rate and strength of adhesion
- Microorganisms can form their own «sticky matrix» for highest persistence and protection = biofilm


Bacteria on stainless steel, scanning electron micrograph of “negative imprint” .

Kalab, M. *et al.* [2008]. *infocus Magazine* 10.

DOI: 10.22443/rms.inf.1.33

# Why is the choice of the sampling device [and the moistening solution] important?


- Force that can be applied for removal of attached cells
- Possibility for sampling of hard-to-reach places [edges, grooves, valves, tubes, crevices, holes etc.]



- Surface area covered by device
- Effectiveness of and capacity for pick-up of cells
- Risk of cross-contamination during sampling
- Survival and/or growth of microorganisms on device during transport
- Release of cells into resuspension buffer/enrichment medium

# Why is the sampling procedure important?


- Optimal moistening solution -> best recovery rate
- Sampled spot/area should be precisely defined  
-> best reproducibility in consecutive samplings [trend analysis]
- **Sampled area/spot should be covered multiple times** with one device  
[going back and forth in two to three different directions]



- **Sampling device should be turned** during sampling to take advantage of its full capacity
- **Sufficient pressure** should be applied without breaking the device  
-> optimal pick-up of adherent cells and biofilm

# Sampling devices available today

## Swabs



## Sponge sticks



## Sponges and wipes



# How do different sampling devices perform?



Recovery rate [% cfu recovered of cfu inoculated]:

Gomez et al. [2012]. J. Food Protection 75: 1077-1082


Device	<i>L. monocytogenes</i> strain 935 (human origin)	<i>L. monocytogenes</i> strain 437/07 (food origin)
Sponge (4 x 8 cm, with 10 ml BPW)	2.81 ± 0.63 BCDEF	0.97 ± 0.36 BCDEF
Premoistened towel (20 x 20 cm)	0.28 ± 0.18 G	0.06 ± 0.03 F
Cotton swab	1.01 ± 0.47 EFG	0.09 ± 0.07 F
Alginate swab	0.22 ± 0.15 G	0.14 ± 0.11 F
Metallic swab	0.24 ± 0.29 G	0.03 ± 0.03 F
Cotton disk	2.53 ± 0.96 CDEFG	1.22 ± 0.67 BCDE
Cotton g Miniroller 1 White microfiber	0.40 ± 0.15 G	0.35 ± 0.20 DEF
MR1S	3.53 ± 1.17 BCD	1.30 ± 0.72 BCD
MR2S Miniroller 2 100% Wool fiber-velour	6.27 ± 1.62 A	1.52 ± 0.37 B
MR3S	3.31 ± 0.10 BCDE	1.41 ± 0.55 BC
MR4S Miniroller 3 100% White polyamide fiber	2.32 ± 0.48 CDEFG	0.43 ± 0.17 CDEF
MR5S	3.69 ± 0.64 BCD	0.59 ± 0.24 BCDEF
MR1V Miniroller 4 White high-density foam	1.79 ± 0.36 DEFG	0.31 ± 0.12 DEF
MR2V	5.05 ± 2.19 AB	1.49 ± 0.65 B
MR3V Miniroller 5 High-density foam flocked with polyamide fiber of 3 mm	4.23 ± 1.53 ABC	0.21 ± 0.18 EF
MR4V	1.90 ± 0.58 CDEFG	0.61 ± 0.33 BCDEF
MR5V	1.33 ± 0.66 DEFG	0.61 ± 0.46 BCDEF
Steel wool pad	0.75 ± 0.44 FG	0.31 ± 0.15 DEF
Petrifilm	0.25 ± 0.29 G	1.20 ± 0.36 BCDE
RODAC ALOA	2.50 ± 1.29 CDEFG	4.15 ± 1.03 A

Used stainless steel table, 10 x 10 cm squares (100 cm<sup>2</sup>), *L. mono.* diluted in peptone water/minced meat suspension, 0.1 ml of 10<sup>6</sup> CFU suspension spread over square, dried for 15 min.

# A new approach to sampling of large surfaces: The N-Light™ MaxiSampler




**The first large surface sampling device  
compatible with rapid, on-site tube tests  
Patent pending**



Firmly attached, robust handle for forceful swabbing,  
removable by simple click mechanism

Giant swab barrel with fine tip, flocked fiber coating (similar to  
standard NEMIS swabs) on outside, hollow inside for dropping of  
N-Light™ tablets to enrichment broth at bottom of N-Light™ test  
tubes





# How to use the N-Light™ MaxiSampler



## 1. PREPARATION

Repeat for 5-6 samples per reservoir and buffer tube



## 2. SAMPLING



Add the N-Light™ buffer to the sterile reservoir.



Mark/label the pouch according to your test plan. Open the sealed pouch.



Remove the N-Light™ Maxi-Sampler from the pouch and wet it in the buffer.



Store the wetted MaxiSampler again in its pouch and use it within 4 hours.



Sample surface areas of up to 0.3 m<sup>2</sup> by scraping the area in 3 directions with some force. Place the N-Light™ Maxi-Sampler back into the pouch after sampling and proceed with step 3 within 4 hours.

### Tips:

- For sampling wet surfaces, using dry MaxiSamplers is acceptable.
- Use the tip of the MaxiSampler swab to sample grooves and holes in the sampling area.
- Change the inclination of the handle during swabbing.

**See also Instructions for Use (inlay in each box of MaxiSamplers, available as PDF)**

# How to use the N-Light™ MaxiSampler



## 3. TRANSFER Procedure depends on the N-Light™ test

<p>N-Light™ <i>Listeria monocytogenes</i> N-Light™ <i>E. coli</i></p>	<p>REF 00009 N-LIGHT™ LISTERIA MONOCYTOGENES</p>	<p>REF 00029 N-LIGHT™ E. COLI</p>	
<p>Open the N-Light™ test tube. Insert the MaxiSampler swab and unclip the swab by turning the handle upwards.</p>	<p>Put the cap back on the test tube and lock the tube [press firmly until you hear a "click"].</p>	<p>Rinsing: Hold tube between thumb and forefinger in horizontal orientation with fiber coating facing down. then incline to one side, then to the other, repeat 10 times</p>	<p>The MaxiSampler swab remains in the tube for the subsequent steps: 4. Incubate 5. Activate 6. Measure</p>

***L. mono., E. coli:*  
MaxiSampler IN**

Important:  
Rinse well with coating facing down before incubation

<p>N-Light™ <i>Salmonella Risk</i> N-Light™ <i>Listeria spp.</i></p>	<p>REF 00014 N-LIGHT™ SALMONELLA RISK</p>	<p>REF 00048 N-LIGHT™ LISTERIA spp.</p>	
<p>Release the MaxiSampler swab into the N-Light™ test tube and carefully screw the cap back onto the tube. DO NOT LOCK THE TUBE.</p>	<p>Rinsing: Hold tube between thumb and forefinger in horizontal orientation with fiber coating facing down. then incline to one side, then to the other, repeat 10 times</p>	<p>Open the test tube again and insert the back of the handle into the MaxiSampler swab. REMOVE THE SWAB FROM THE TUBE. Discard the swab.</p>	<p><i>Salmonella Risk only:</i> Add one starter tablet to the test tube and briefly shake the tube. All: NOW LOCK THE TUBE.</p>

***Salmonella Risk, L. spp.:*  
MaxiSampler OUT**

Important:  
Rinse well with coating facing down before removal and incubation

(See also IFU)

# Recovery rate of N-Light™ MaxiSampler

- Bacterial cells diluted in sterile broth, 10x 10 µl «spotted» on 121 cm<sup>2</sup> stainless steel plates
- 1.5 h drying at room temperature
- MaxiSamplers moistened with N-Light™ Neutralizer, scrub samplers used as received.
- Surface swabbed in 3 directions while turning devices
- Devices resuspended in equal volumes of buffered peptone water (BPW)
- Inoculated and recovered (resuspended) cells [colony forming units] counted on plate count agar

Sampling Device	Test strain	CFU inoculated	CFU recovered after sampling	Recovery rate
N-Light™ MaxiSampler	<i>Listeria monocytogenes</i> ATCC 19111	39'000 ± 3'606	5'262 ± 127	14% ± 1%
	<i>Salmonella</i> Typhimurium ATCC 14028	17'000 ± 3'754	5'201 ± 585	30% ± 8%
3M Scrub Sampler	<i>Listeria monocytogenes</i> ATCC 19111	39'000 ± 3'606	6'319 ± 2'237	16% ± 6%
	<i>Salmonella</i> Typhimurium ATCC 14028	17'000 ± 3'754	3'606 ± 321	21% ± 5%

# Use of N-Light™ MaxiSampler with the N-Light™ *Listeria monocytogenes* Test



- Bacterial cells diluted in BHI broth, 10x 10 µl «spotted» on stainless steel (121 cm<sup>2</sup>)
- 1 h drying at room temperature
- MaxiSamplers moistened with N-Light™ Neutralizer, stainless steel plates swabbed in 3 directions
- MaxiSamplers left IN, incubation and read-out according to IFU for N-Light™ *Listeria monocytogenes*

Inoculation level [cfu/ 100 cm <sup>2</sup> steel plate]	Positive N-Light™ <i>Listeria monocytogenes</i> tests	
	<i>L. mono.</i> ATCC 13932	<i>L. mono.</i> ATCC 19111
1000-2000	5 / 5	5 / 5
100-200	5 / 5	5 / 5
10-20	4 / 5	5 / 5

# N-Light™ MaxiSampler compared to N-Light™ swabs BP

N-Light *L. mono.* tests

Moistening solution:  
N-Light™ Neutralizer

MaxiSampler/swabs  
left IN

Threshold for positive  
tests: 10'000 RLU

<b>MaxiSampler</b>	<b>RLU</b>	<b>RLU</b>	<b>Positives</b>
Stainless steel, 11 x 11 cm	average	st. dev.	
<i>L. mono.</i> ATCC 13932, 12'000 cfu	391'263	81'293	<b>5/5</b>
<i>L. mono.</i> ATCC 13932, 1'200 cfu	612'271	329'705	<b>5/5</b>
<i>L. mono.</i> ATCC 13932, 120 cfu	252'008	237'596	<b>4/5</b>
<i>L. mono.</i> ATCC 13932, 12 cfu	2'775	1'167	<b>0/5</b>
Negative control (Sterile BHI on steel)	1'878	NA	<b>0/1</b>
<b>Swabs BP</b>	<b>RLU</b>	<b>RLU</b>	<b>Positives</b>
Stainless steel, 2.5 x 2.5 cm	average	st. dev.	
<i>L. mono.</i> ATCC 13932, 12'000 cfu	3'026'692	479'233	<b>5/5</b>
<i>L. mono.</i> ATCC 13932, 1'200 cfu	3'151'930	962'436	<b>5/5</b>
<i>L. mono.</i> ATCC 13932, 120 cfu	3'029'002	1'152'811	<b>5/5</b>
<i>L. mono.</i> ATCC 13932, 12 cfu	2'100'827	1'922'773	<b>3/5</b>
Negative control (Sterile BHI on steel)	2'759	NA	<b>0/1</b>

-> N-Light™ swabs BP perform better than MaxiSamplers when used for small areas ( $\leq 10 \text{ cm}^2$ )

# Use of N-Light™ MaxiSampler with the N-Light™ *Escherichia coli* Test




- Bacterial cells diluted in TSB, 10x 10 µl «spotted»
- 2 h drying
- MaxiSamplers with cells kept 2 h at room temp.
- Incubation and read-out according to IFU for N-Light™ *E. coli*

Negative control:  
 (sterile broth on steel)  
 MaxiS IN 1'008 RLU  
 MaxiS OUT 1'465 RLU

Threshold for positives: 10'000 RLU

Stainless steel, 25-30 cfu inoculated per 121 cm<sup>2</sup> square



# Use of N-Light™ MaxiSampler with the N-Light™ *Salmonella* Risk Test

Similar experimental procedure as for the N-Light *E. coli* test, except for removal of MaxiSampler after rinsing, moistening with Neutralizer



Inoculation level <i>Salmonella</i> Typhimurium ATCC 14028 [cfu/ 121 cm <sup>2</sup> steel plate]	Result N-Light™ <i>Salmonella</i> Risk Test	
	Positives	RLU, average ± SD
790	3 / 3	538'738 ± 202'228
70	3 / 3	206'659 ± 38'642
0 [sterile TSB broth]	0 / 2	1'106

Threshold for positives: 20'000 RLU

# Use of N-Light™ MaxiSampler with the N-Light™ *Listeria* spp. Test



Cells diluted in BHI broth, otherwise similar experimental procedure as for the *E. coli* and *Salmonella* Risk tests, MaxiSamplers removed after rinsing [before incubation]

Inoculation level and strain [cfu/ 121 cm <sup>2</sup> steel plate]	Result N-Light™ <i>Listeria</i> spp. tests	
	Positives	RLU
98 cfu <i>L. mono.</i> ATCC 13923	2 / 2	11'173 / 15'534
230 cfu <i>L. innocua</i> ATCC 13923	2 / 2	24'097 / 17'703
0 [sterile BHI]	0 / 5	1'147 / 1'094

Threshold for positives: 4'000 RLU



# Field tests with N-Light™ MaxiSampler

Ready-to-eat food factory 1

Moistening solution	Positives N-Light™ <i>Salmonella</i> Risk	PCR analysis of N-Light tests, <i>Salmonella</i> positive
N-Light™ BPW	4 / 20 [20%]	0 / 20
N-Light™ Neutralizer	5 / 20 [25%]	0 / 20

Ready-to-eat food factory 2

Positives N-Light™ <i>L. mono.</i>	ISO analysis of developed N-Light tests, <i>L. mono.</i> positive
2 / 23 [9%]	0 / 23



# Summary

- **Sampling** is determines to a large part the **quality of environmental monitoring data**
- **Microorganisms** can be **firmly attached to surfaces** (and hard to pick-up), particularly in **biofilms**
- There is **no one-size fits all sampling device**
- Large [MaxiSampler, sponges, sponge sticks] and small [swabs] sampling devices each have **distint advantages and disantvantages**
- The **surface contact material of sampling devices matters**, but more research is needed
- The novel **NEMIS N-Light™ MaxiSampler** allows **efficient sampling of large surfaces**
- The MaxiSampler is **compatible with all N-Light™ pathogen/indicator tests**
- **Overloading of MaxiSamplers with food matrix/dirt should be avoided**



Thanks for your  
participation !